

The Purple Stirrup™

The Innovative KYRA COMFORT™ 500 Stirrup



The KYRA COMFORT™ 500 Stirrup The KYRA COMFORT™ 500 SB Stirrup

- **BETTER INFECTION CONTROL:** Optional KYRA SECURE-LOK™ Closure Straps¹ eliminate VELCRO® brand fastener straps and are designed to improve infection control.
- **UNPARALLELED PELVIC SITE ACCESS FOR ROBOTIC ASSISTED PROCEDURES:** Wide lithotomy range of motion from +85° to -55°. The unique patent pending design provides 22° more lithotomy range of motion than traditional stirrups.



OPTIONAL SECURE-LOK™
CLOSURE STRAPS



140° LITHOTOMY RANGE OF MOTION



innovative solutions for today's modern OR™

The KYRA COMFORT™ 500 Stirrup

- **NO MORE LOST CLAMPS:** The Integrated KYRA SECURE-RELEASE™ Clamp¹ solves the “lost clamp problem”. If the clamp is damaged it can be easily replaced or removed between cases allowing uninterrupted use.
- **ACCOMMODATES VIRTUALLY ALL PATIENTS:** The larger/wider KYRA COMFORT™ 500 Stirrup boots and the standard sized boot on the KYRA COMFORT™ SB Stirrup accommodate virtually all patient leg sizes.
- **MORE PRESSURE REDUCING PADDING:** Kyra Medical’s unique stirrup design¹ allows ½” more padding for pressure relief than traditional stirrups². Unlike traditional boot pads, Kyra’s pad protects the medial and lateral aspect of the ankle.
- **ELIMINATES DANGEROUS PINCH HAZARDS³:** KYRA COMFORT Stirrup’s curved support rods¹ eliminate dangerous pinch hazards that traditional stirrups present.



FLEXIBLE BOOT ACCOMMODATES DIFFERENT LEG SIZES



NO MORE LOST CLAMPS! OPTIONAL INTEGRATED CLAMP

PRODUCT:	VELCRO® Brand Fastener Straps*	SECURE-LOK™ Closure Straps	SECURE-RELEASE™ Integrated Clamps
KYRA COMFORT™ 500 Stirrups w/ traditional straps PC: KYRA5000	YES	NO	NO
KYRA COMFORT™ 500 Stirrups w/SECURE-LOK™ straps PC: KYRA5100	NO	YES	NO
KYRA COMFORT™ 500 Stirrups w/ integrated clamps & traditional straps PC: KYRA5500	YES	NO	YES
KYRA COMFORT™ 500 Stirrups w/ integrated clamps & SECURE-LOK™ straps PC: KYRA5600	NO	YES	YES
KYRA COMFORT™ 500 SB Stirrups w/ SECURE-LOK™ straps PC: KYRA5105	NO	YES	NO
KYRA COMFORT™ 500 SB Stirrups w/ integrated clamps & SECURE-LOK™ straps PC: KYRA5605	NO	YES	YES


Specifications

The KYRA COMFORT™ 500 boot accommodates up to the 95th percentile leg size patients weighing 500 lb.

The KYRA COMFORT™ 500 SB boot accommodates heavier patients with smaller legs (patients weighing up to 500 lb.).

Together the KYRA COMFORT™ 500 and the KYRA COMFORT™ 500 SB accommodate 95% of patients weighing 100 -up 500 lbs.

Patient Weight Capacity: 500 lbs (227 kg) Abduction Range: +25° to -9°
Lithotomy Range: +85° to -55°



THE PURPLE STIRRUP™ CART

The Purple Stirrup™ Cart allows for convenient and safe storage of any stirrup.



innovative solutions for today's modern OR™

888-611-KYRA (5972)

kyramedical.com