



KYRA[®] 800 STIRRUP

STRONG. LIGHTWEIGHT. LARGE BOOT

The KYRA Comfort™ 800 stirrup is designed for patients with a BMI of greater than 29 with a large flexible lateral fin and can support patients up to an 800 lbs/363 kg weight capacity.

KYRA boots will safely accept larger legs with 1" of supportive boot padding.

All Comfort 800's come standard with the durable Secure-Lok™ straps which help reduce infection control and the patented Core-Grip™ technology which makes the stirrup stronger and lighter vs other stirrup brands.



Part Numbers

KYRA8100 KYRA 800 Stirrups w/ Secure-Lok™ straps

Specifications

- KYRA Comfort 800 weight: 14.6 lbs/6.6 kg each
- More working space with wider abduction range and boot rotational angle
- Lithotomy Range: +85° to -55°, Abduction Range: +25° to -9°, Boot Rotational Angle: +/-25°
- Lightweight and strong with KYRA's Core-Grip™ technology

