



# KYRA<sup>®</sup> 250 STIRRUP

## NON-FINNEO. JUNIOR.

The KYRA 250 is a junior stirrup designed for patients aged 6-12 years old but also effective in smaller/lighter adults weighing up to 250 lbs/113 kg.

The non-finned boot has ¾" of supportive boot padding and larger ankle well for more secure fit.

All KYRA 250's come standard with the durable Secure-Lok™ straps which help reduce infection control and the patented CoreGrip™ technology which makes the stirrup stronger and lighter vs other stirrup brands.



### Part Numbers

KYRA7180 KYRA 250 Stirrups w/ Secure-Lok™ straps

### Specifications

- KYRA 250 stirrup weight: 10.2 lbs/4.6 kg each
- More working space with wider boot rotational range
- Lithotomy Range: +85° to -55°, Abduction Range: +25° to -9°, Boot Rotational Angle: +/-25°
- Lightweight and strong with KYRA's CoreGrip™ technology

