



KYRA[®] 350 STIRRUP

STRONG. LIGHTWEIGHT.

The KYRA Comfort™ 350 stirrup is designed for patients with a BMI of less than 29 with a large flexible lateral fin and can support patients up to 350 lbs/159 kg weight capacity.

KYRA boots will safely accept larger legs with $\frac{3}{4}$ " of supportive boot padding.

All Comfort 350's come standard with the durable Secure-Lok™ straps which help reduce infection control and the patented CoreGrip™ technology which makes the stirrup stronger and lighter vs other stirrup brands.



Part Numbers

KYRA3100 KYRA 350 Stirrups w/ Secure-Lok™ straps

Specifications

- KYRA Comfort 350 weight: 13.2 lbs/6 kg each
- More working space with wider abduction range and boot rotational angle
- Lithotomy Range: +85° to -55°, Abduction Range: +25° to -9°, Boot Rotational Angle: +/-25°
- Lightweight and strong with KYRA's CoreGrip™ technology

