



# KYRA<sup>®</sup> 400 STIRRUP

## NON-FINNEO. ROBOT FRIENDLY.

The KYRA 400 or "KYRA Blue" is a stirrup designed for customers that prefer a non-finned boot. KYRA 400 stirrups can support patients up to a 400 lbs/181 kg weight capacity.

It has a premium stainless-steel gas shock and a boot with 3/4" of supportive boot padding.

All KYRA 400's come standard with the durable Secure-Lok<sup>™</sup> straps which help reduce infection control and the patented CoreGrip<sup>™</sup> technology which makes the stirrup stronger and lighter vs other stirrup brands.



### Part Numbers

KYRA7100 KYRA 400 Stirrups w/ Secure-Lok<sup>™</sup> straps

### Specifications

- KYRA 400 weight: 12.2 lbs/5.5 kg each
- More working space with wider abduction and boot rotational range
- Lithotomy Range: +85° to -55°, Abduction Range: +25° to -9°, Boot Rotational Angle: +/-25°
- Lightweight and strong with KYRA's CoreGrip<sup>™</sup> technology

