

KYRA[®] Spine

Radiolucent. Fabric Free Straps. 600 lbs. Capacity.



- INTEGRATED GEL PADS
- ALLOWS FOR MORE ABDOMEN SPACE



- DURABLE SILICONE STRAPS
- HOOK & LOOP FREE

KYRA Curve[™] Spine Frame System KYRA9000K

Provides a stable platform for spine surgeries including laminectomies, decompressions and disc surgery. Intraoperatively adjustable to the desired height to get optimal surgical site exposure. C-arm compatible. Fits spine table rails and tabletops. Unique silicone straps provides better cleaning vs. hook & loop straps. Pad set included.

Patient Weight Capacity: 600 lbs. (272 kg)

Device Weight: 24.4 lbs. (11 kg)

WHAT'S INCLUDED



KYRA9000
Curve Spine Frame
1 Included

+



KYRA9010
Curve Crank Handle
1 Included

+



KYRA1100
Curve Spine Frame Cart
1 Included

+




KYRA9030
Secure-Lok[™] Strap Connector
For use with tables without
side rails.

WHAT'S OPTIONAL



KYRA9020
Curve Spine Frame
Comfort Covers Set

KYRA® Curve™ Spine Frame Disposables



KYRA Curve Comfort Cover Set KYRA9020 (6 pack)

Disposable pack for use with KYRA Curve Spine Frame. Includes cosmetic foam pad covers, head block, lateral rod pad and arm cushions.



KYRA Curve Standard Cover Set KYRA9021 (6 pack)

Disposable pack for use with KYRA Curve Spine Frame. Includes standard foam pad covers, head block, lateral rod pad and arm cushions.



KYRA Foam Head Block KYRA9022 (6 pack)

Provides support, stability, elevation and safe alignment of the head and neck.