




The KYRA® Comfort™ Family of Stirrups

Maximum range of motion for precision positioning.

Unique curved rod design engineered to support the widest range of lithotomy and abduction—ideal for optimizing positioning, surgical site access and robotics.

- ✓ Three weight capacity options to accommodate patients up to 800 lbs (363 kg).
- ✓ Easy-clean, durable, fabric-free silicone straps for improved infection control.
- ✓ Larger boot design with pressure-reducing padding accommodates a wider range of patient anatomies.
- ✓ Engineered for maximum up-time with on-site repairs with KYRAComplete™ Repair Kits.




Learn more about the **Comfort™ Family** →


Learn more at kyramedical.com
Your trusted partners in patient positioning



The KYRA Difference




✓ Curved rod provides larger range of motion optimized for low lithotomy and robotics



✓ Free-floating boot design safely maintains proper angle for patient's leg



✓ Flexible lateral fin, increased padding and larger boot for peroneal nerve protection



✓ Easy-clean, fabric-free silicone SecureLok™ straps help improve infection control and snap in securely with audible click



✓ Additional 5° of boot rotation internally and externally for proper alignment on wider tables

✓ Maximize uptime with repair kits, extended warranty and premium servicing

Find Your Fit

Model Comparison Guide



	Comfort 350 Stirrups	Comfort 500 Stirrups	Comfort 800 Stirrups
PATIENT WEIGHT CAPACITY	350 lbs. (159 kg)	500 lbs. (227 kg)	800 lbs. (363 kg)
SUGGESTED PATIENT BMI	Under 29	29 or greater	29 or greater
STIRRUP WEIGHT	13.2 lbs. (6 kg)	14.6 lbs. (6.6 kg)	14.6 lbs. (6.6 kg)
BOOT PADS THICKNESS	¾" (19 mm) pads included	1" (25 mm) pads included	1" (25 mm) pads included
LITHOTOMY RANGE	+85° to -55° (22° more than competitive stirrups)		
USE CASE	Optimized for low lithotomy and robotics cases	Optimized for low lithotomy and robotics cases	Optimized for low lithotomy and robotics cases with bariatric patients
PRODUCT NO.	KYRA3100 - Silicone Straps	KYRA5100 - Silicone Straps	KYRA8100 - Silicone Straps

All Comfort Stirrups are sold in pairs with pads included. Requires two Secure Release™ Clamps.

Products may be covered by one or more issued or pending patents. Visit kyramedical.com/innovation for details.

Scan to learn more at kyramedical.com
Your trusted partners in patient positioning

