

NEW!



KYRA Connect™

SPINE SYSTEM

A space-saving mobile unit that transforms your OR table into an operative spine & imaging platform

Bringing unmatched functionality and flexibility to the OR table.



Explore the KYRA Connect™ Spine System →

KYRA Connect™ SPINE SYSTEM

Ideal for overflow cases and increasing capacity.

The Kyra Connect Spine System transforms most standard OR table into a versatile spine and imaging platform—without the cost or footprint of a traditional spine system.

- ✓ **Compact Functionality** — “Jackson-style” table performance with the footprint of an IV-pole.
- ✓ **Imaging Compatible**—Full O-arm and C-arm access.
- ✓ **Modular Design**—Configurable to support procedural requirements.
- ✓ **Designed for Clinicians**—Engineered for simple setup and storage.



Specifications:
Patient Weight Capacity: 500 lbs. (227 kg)
Product Weight: 100 lbs. (50 kg)

Products covered by one or more issued or pending patents.

Discover our multi-patented design:



Table Mounts securely attach Connect rails to existing OR table



Patented 5-axis flex joint supports range of movement of the table without restriction



Integrated bubble level to assist in leveling table



Powered column provides optimal surgical and imaging height

Lightweight, mobile system with footprint of an IV pole easily transports and stores accessories



Auto-locking casters for ease of use and enhanced security



Simple handheld remote controls locking, storage and height adjustments

Choose the Components. Control the Setup.

Multiple configurations support a wide range of spine procedures while maintaining imaging access. See the back page for a recommended setup.



Compatible with the Curve™ Spine Frame



Supports prone, lateral & supine positioning.

A Recommended Setup, Adaptable by Design

Mix and match accessories to create the exact configuration your spine procedures require. Talk with a specialist to build your ideal Kyra Connect configuration.



A **Connect™ Spine System**
KYRA9500

B **Prone Head Support System**
KYRA9071K
Prone Face Cushion, Adult
KYRA9080

C **Connect™ Head Plate**
KYRA9050

D **Connect™ Swivel Armboard**
KYRA2030

E **Connect™ Rail Clamp, 6"**
KYRA9525

F **Connect™ Prone Support Set**
KYRA9052K

Connect™ Comfort Covers
KYRA9063

G **Connect™ Table Mount Clamps**
KYRA9511 (US)

H **Connect™ Supine Top, Kit of 4**
KYRA9056K

I **Connect™ Power Cord (not shown)**
KYRA9540 (US)

Connect with us to customize your ideal Kyra Connect configuration.