

# KYRA Curve™ Spine Frame System

## Optimized positioning for spine procedures

The KYRA Curve Spine Frame provides radiolucent, optimal positioning across a wide range of spine procedures. Adjustable height and width enhances intervertebral access and reduces abdominal pressure for a range of patient anatomies.

- ✓ Ideal for spine procedures requiring reverse lordosis including laminectomies, decompressions, fusions and discectomies.
- ✓ Intraoperative height adjustment for optimal surgical site exposure.
- ✓ C-arm compatible and fits dual-rail spine systems, standard OR tables and imaging tables.
- ✓ Laterally adjustable pads allow abdomen to hang freely.
- ✓ Easy-to-clean silicone straps are durable while improving infection control.
- ✓ Two size options ensure neutral alignment of the head and neck for full range of patient heights.



Crank handle allows for adjustable intraoperative positioning



Learn more about the Curve & Mini Curve →

Learn more at [kyramedical.com](http://kyramedical.com)  
Your trusted partners in patient positioning





# Curve™ Spine Frame System

KYRA9000K - Standard  
KYRA9005K - Mini

Pad set included.

**Patient Weight Capacity:** Standard - 600 lbs. (272 kg) • Mini - 500 lbs. (227 kg)

**Device Weight:** Standard - 24.4 lbs. (11 kg) • Mini - 22 lbs. (10 kg)

Fits “Jackson-style” rails, standard OR and imaging tables

## One Design. Two Sizes. Optimal Positioning.

The only ‘Wilson-style’ frame on the market designed for shorter patients; the KYRA Mini Curve™ Spine Frame offers all the benefits of the KYRA Curve™ Spine Frame with a 30% shorter frame to accommodate patients of 5’9” and under. See sizing chart below.

### WHAT'S INCLUDED

Each system includes one Spine Frame (Standard or Mini) with pads plus Crank Handle, Storage Cart and Curve Strap Buckles.




#### Standard Curve Spine Frame

For patients 5’7” and above, up to 600 lbs. (272 kg).

KYRA9000K

OR



#### Mini Curve Spine Frame

For patients 5’9” and under, up to 500 lbs. (227 kg).

KYRA9005K



#### Curve Crank Handle

One rotational handle for controlled, secure lordosis positioning.

KYRA9010

+



#### Curve Spine Frame Cart

Safely stores and transports spine frame and components.

KYRA1100 (Standard)  
KYRA1110 (Mini)

+



#### Curve Strap Buckle (Pair)

For use with tables without side rails.

KYRA9030

### WHAT'S OPTIONAL


#### Curve™ Comfort Covers

KYRA9020S Standard  
KYRA9024 Mini

6 sets per pack



- ✓ Layer of supportive foam absorbs excess moisture, dissipates heat and helps reduce risk of skin shear.
- ✓ Single-use for improved infection control.



\* Using correct size frame ensures neutral alignment of the head and neck during prone positioning.

Products may be covered by one or more issued or pending patents. Visit [kyramedical.com/innovation](http://kyramedical.com/innovation) for details.

Scan to learn more at [kyramedical.com](http://kyramedical.com)

Your trusted partners in patient positioning

