



KYRA® Swivel Armboards

From simple rotation to imaging-ready prone positioning

KYRA swivel armboards provide a variety of adaptable upper extremity supports for precise surgical positioning. With smooth movement and reliable stability, these solutions enable efficient positioning throughout the procedure.

- ✓ Controlled movement for accurate arm placement
- ✓ Multiple options to support a variety of procedural needs and clinical preferences
- ✓ Silicone cover prevents fluid ingress and protects hardware



Learn more at kyramedical.com
Your trusted partners in patient positioning



Standard Swivel Armboards



Simple, table-mounted rotation for efficient positioning

- ✓ Two swivel post sizes for optimized positioning range.
- ✓ Multi-task armboard, typically used for lateral decubitus and beachchair positions.
- ✓ 1" (25 mm) thick Clamshell Pad included for comfort and protection.

Connect™ Swivel Armboard



Low-profile design reduces imaging interference

- ✓ Designed for optimal positioning without C-arm interference.
- ✓ Lightweight, low-profile design.
- ✓ Integrated SocketLok™ clamp for quick table attachment.
- ✓ Clamshell Pad included for comfort and protection.

Find Your Fit

Model Comparison Guide

	🌀 Swivel, Small	🌀 Swivel, Large	🌀 Connect Swivel
POSITIONING RANGE	Universal	Universal	Prone
ARMBOARD PAD STYLE	Clamshell Pad, 1" thick	Clamshell Pad, 1" thick	Clamshell Pad, 1" thick
PAD INCLUDED	Yes	Yes	Yes
POST SIZE	16" H x 6.5"W (40.6 cm x 16.5 cm)	18.2"H x 12.5"W (46.2 cm x 31.7 cm)	n/a
INTEGRATED CLAMP	No	No	Yes
PRODUCT NO.	KYRA2304	KYRA2300	KYRA2302-US
USE CASE	Basic support, e.g. well arm in beachchair	Multi-task flexibility, e.g. lateral positioning	Maximum positioning flexibility with single-handed adjustment

Need quick, single-handed positioning?
Ask about our traditional

Trigger-style Armboards



Scan to learn more at [kyramedical.com](https://www.kyramedical.com)

Your trusted partners in patient positioning